

New England Biolabs Product Specification

Product Name: SARS-CoV-2 LAMP Primer Mix (N/E)
Catalog #: S1883B
Concentration: 10X Concentrate
Shelf Life: 24 months
Storage Temp: -20°C
Composition (1X): Proprietary
Specification Version: PS-S1883B v1.0
Effective Date: 15 Jul 2020

Assay Name/Specification (minimum release criteria)

Endonuclease Activity (Nicking) - A 50 µl reaction in NEBuffer 2 containing 1 µg of supercoiled PhiX174 DNA and a minimum of 5 µl of SARS-CoV-2 LAMP Primer Mix (N/E) incubated for 4 hours at 37°C results in <10% conversion to the nicked form as determined by agarose gel electrophoresis.

Non-Specific DNase Activity (16 Hour) - A 50 µl reaction in NEBuffer 2 containing 1 µg of T3 or T7 DNA in addition to a reaction containing Lambda-HindIII DNA and a minimum of 5 µl of SARS-CoV-2 LAMP Primer Mix (N/E) incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

Phosphatase Activity (pNPP) - A 200 µl reaction in 1M Diethanolamine, pH 9.8, 0.5 mM MgCl₂ containing 2.5 mM *p*-Nitrophenyl Phosphate (pNPP) and a minimum of 20 µl of SARS-CoV-2 LAMP Primer Mix (N/E) incubated for 4 hours at 37°C yields <0.0001 unit of alkaline phosphatase activity as determined by spectrophotometric analysis.

RNase Activity (Extended Digestion) - A 10 µl reaction in NEBuffer 4 containing 40 ng of a 300 base single-stranded RNA and a minimum of 1 µl of SARS-CoV-2 LAMP Primer Mix (N/E) is incubated at 37°C. After incubation for 16 hours, >90% of the substrate RNA remains intact as determined by gel electrophoresis using fluorescent detection.

*"One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information."*



Date 15 Jul 2020

Derek Robinson
Director of Quality Control

