

New England Biolabs Product Specification

Product Name: CLIP-Cell™ Block
Catalog #: S9220S
Shelf Life: 60 months
Storage Temp: -20°C
Specification Version: PS-S9220S v3.0
Effective Date: 11 Apr 2019

Assay Name/Specification (minimum release criteria)

Identity (Mass Spectrometry) - The observed molecular mass of CLIP-Cell™ Block is 285.0 Da +/- 1 Da as determined by mass spectrometry analysis.

In Vitro Protein Labeling (Blocking Assay) - A 50 µl reaction in 1X PBS and 1 mM DTT containing 5 µM of CLIP-tag® Purified Protein and a minimum of 20 µM of CLIP-Cell™ Block for 30 minutes followed by 10 µM CLIP-Cell™ TMR-Star is incubated for 1 hour at 37°C results in the expected absence of labeled product that is visualized on SDS-PAGE by fluorescent detection.

Physical Purity (HPLC) - CLIP-Cell™ Block is ≥ 90% pure as determined by HPLC analysis.



Date 11 Apr 2019

Derek Robinson
Director of Quality Control

