

Citations for IMPACT

Protein Purification

2010

Collins, E. D., et. al. (2010) Cloning the human vitamin D receptor into the pTwin-1 expression vector. *J Steroid Biochem Mol Biol* (pTWIN1) PubMed ID: [20171280](#)

Lu, M., et. al. (2010) Purification of untagged HIV-1 reverse transcriptase by affinity chromatography. *Protein Expr. Purif.* 71, 231-239.

Manconi, B., et. al. (2010) Expression, purification, phosphorylation and characterization of recombinant human statherin. *Protein Expr Purif.* 69, 219-225.

Yu, H. H., et. al. (2010) Expressed protein ligation for the preparation of fusion proteins with cell penetrating peptides for endotoxin removal and intracellular delivery. *Biochim Biophys Acta*, (IN PRESS) (pTXB1) PubMed ID: [20170629](#)

2009

Che, N., et. al. (2009) Soluble expression and one-step purification of a neurotoxin Huwentoxin-I in *Escherichia coli*. *Protein Expr Purif.* 65, 154-159.

2008

Babu, S.K., et. al. (2008) Construction of intein-mediated hGMCSF expression vector and its purification in *Pichia pastoris*. *Protein Exp. Purif.* 57, 201-205.

Chen, Y. Q., et. al. (2008) Expression of a cytotoxic cationic antibacterial peptide in *Escherichia coli* using two fusion partners. *Protein Exp. Purif.* 57, 303-311.

Zhao, Z., et. al. (2008) Purification of green fluorescent protein using a two-intein system. *Appl. Microbiol. Biotechnol.* 77, 1175-1180.

2007

Bhat, R. K. and Berger, S. (2007) New and easy strategy for cloning, expression, purification, and characterization of the 5S subunit of transcarboxylase from *Propionibacterium f. shermanii*. *Prep. Biochem. Biotechnol.* 37, 13-26.

Diao, H., et. al. (2007) Intein-mediated expression is an effective approach in the study of beta-defensins. *Biochem Biophys Res Commun.* 357, 840-846.

Hassiepen, U., et. al. (2007) A sensitive fluorescence intensity assay for deubiquitinating proteases using ubiquitinrhodamine110-glycine as substrate. *Anal Biochem.* 371, 201-207.

Stasser, J. P., et. al. (2007) A multinuclear copper(I) cluster forms the dimerization interface in copper-loaded human copper chaperone for superoxide dismutase. *Biochem.* 46, 11845-11856.

2006

Demidov, V.V., and Broude, N.E. (2006) Profluorescent protein fragments for fast bimolecular fluorescence complementation *in vitro*. *Nat. Protoc.* 1, 714-719.

Esipov, R.S., et. al. (2006) Production and purification of recombinant human glucagons overexpressed as intein fusion protein in *Escherichia coli*. *Protein Pept. Lett.* 13, 343-347.

Hackenberger, C. P., et. al. (2006) Expression of N-terminal Cys-protein fragments using an intein refolding strategy. *Bioorg. Med. Chem.* 14, 5043-5048.

Wojtaszek, J., et al. (2006) LT CI, a novel chymotrypsin inhibitor of the potato I family from the earthworm *Lumbricus terrestris*. Purification, cDNA cloning, and expression. *Comp. Biochem. Physiol B. Biochem. Mol. Biol.* 143, 465-472. (pTYB11)

2005

Ge, X., et. al. (2005) Self-cleavable stimulus responsive tags for protein purification without chromatography. *J. Am. Chem. Soc.* 127, 11228-9.

Goodin, J. L., et al. (2005) *Yersinia pestis* outer membrane type III secretion protein YscC: expression, purification, characterization, and induction of specific antiserum. *Protein Expr. Purif.* 40, 152-163. (pTYB11)

Ingham, A.B., et al. (2005) A versatile system for the expression of nonmodified bacteriocins in *Escherichia coli*. *J. Appl. Microbiol.* 98, 676-683. (pTYB1)

Morassutti, C., et al. (2005) Expression of SMAP-29 cathelicidin-like peptide in bacterial cells by intein-mediated system. *Protein Expr. Purif.* 39, 160-168. (pCYB1, pTYB11)

Rebets, Y., et al. (2005) DNA -binding activity of LndI protein and temporal expression of the gene that upregulates landomycin E production in *Streptomyces globisporus* 1912. *Microbiology*, 281-290. (pTYB2)

2004

Grzybowska, B., Szweda, P. and Synowiecki, J. (2004) Cloning of the thermostable alpha-amylase gene from *Pyrococcus woesei* in *Escherichia coli*: isolation and some properties of the enzyme. *Mol. Biotechnol.* 26,101-110. (pTYB2)

Guo, C., et al. (2004) Intein-mediated fusion expression, high efficient refolding, and one-step purification of gelonin toxin. *Protein Expr Purif.* 37, 361-367.

2003

Bonsager, B.C., et al. (2003) Purification and characterization of the beta-trefoil fold protein barley alpha-amylase/subtilisin inhibitor overexpressed in *Escherichia coli*. *Protein Expr. Purif.* 2, 185-93. (pTYB1)

Salazar, J.C., et al. (2003) Coevolution of an aminoacyl-tRNA synthetase with its tRNA substrates. *Proc. Natl. Acad. Sci. USA.* 100, 13863-13868. (pTYB1)

Schaffer, L., et al. (2003) Assembly of high-affinity insulin receptor agonists and antagonists from peptide building blocks. *Proc. Natl. Acad. Sci. USA.* 100, 4435-9.

Wuebbens, M. M. and Rajagopalan, K. V. (2003) Mechanistic and mutational studies of *Escherichia coli* molybdopterin synthase clarify the final step of molybdopterin biosynthesis. *J. Biol. Chem.* 278, 14523-14532. (pTYB3)

2002

Humphries, H. E., Christodoulides, M. and Heckels, J. E. (2002) Expression of the class 1 outer-membrane protein of *Neisseria meningitidis* in *Escherichia coli* and purification using a self-cleavable affinity tag. *Protein Expr. Purif.* 26, 243-8. (pTWIN 1)

Moolenaar GF, et al. (2002) Cho, a second endonuclease involved in *Escherichia coli* nucleotide excision repair. *Proc. Natl. Acad. Sci. USA.* 99, 1467-1472. (pTYB1)

Nakano, S., et al. (2002) Multiple pathways of Spx (YjbD) proteolysis in *Bacillus subtilis*. *J. Bacteriol*, 184, 3664-3670. (pTYB1, pTYB2)

Singleton, S. F., et al. (2002) Intein-mediated affinity-fusion purification of the *Escherichia coli* RecA protein. *Protein Expr. Purif.* 26, 476-488. (pTXB3)

Wojciak, J.M., et al. (2002) Arm-site binding by lambda -integrase: solution structure and functional characterization of its amino-terminal domain. *Proc. Natl. Acad. Sci. USA*. 99, 3434-3439. (pTYB11)

Wojciechowski CL, Cardia JP, and Kantrowitz ER. (2002) Alkaline phosphatase from the hyperthermophilic bacterium *T. maritima* requires cobalt for activity. *Protein Sci*. 4, 903-911. (pTYB12)

Zhong, Q., et al. (2002) Endosomal localization and function of sorting nexin 1. *Proc. Natl. Acad. Sci. USA*. 99, 6767-6772.

2001

Dorsch, S., et al. (2001) The VP1-unique region of parvovirus B19: amino acid variability and antigenic stability. *J. Gen. Virol*. 82 (pt 1) 191-199.

Hanzawa, H., et al. (2001) *In vitro* assembly of phytochrome B apoprotein with synthetic analogs of the phytochrome chromophore. *Proc. Natl. Acad. Sci. USA*. 98, 3612-3617. (pTYB2)

Protein labeling, ligation and cyclization

2010

Chen, J., et al. (2010) Chemically ubiquitylated PCNA as a probe for eukaryotic translesion DNA synthesis. *Nat. Chem. Biol*. 6, 270-272. (pTYB1) PubMed ID: [20208521](#).

Karukurichi, K. et al. (2010) Analysis of p300/CBP histone acetyltransferase regulation using circular permutation and semisynthesis. *J. Am. Chem. Soc*. 132, 1222- 1223. (pTYB2) PubMed ID: [20063892](#).

Lu, X., et al. (2010) Designed semisynthetic protein inhibitors of ub/ubl e1 activating enzymes. *J. Am. Chem. Soc*. 132, 17481749. (pTXB1) PubMed ID: [20099854](#).

2009

Hondal, R. J. (2009) Using chemical approaches to study selenoproteins-Focus on thioredoxin reductases. *Biochim. Biophys. Acta*. 1790, 1501-1512. PubMed ID: [19406205](#).

Chattopadhyaya, S., et al. (2009) Use of intein-mediated protein ligation strategies for the fabrication of functional protein arrays. *Methods Enzymol*. 462, 195-223. PubMed ID: [19632476](#).

Chattopadhyaya, S., et al. (2009) Site-specific covalent labeling of proteins inside live cells using small molecule probes. *Bioorg. Med. Chem.* 17, 981-9. (pTWIN)
PubMed ID: [18261914](#).

2008

Skrisovska, L. et al. (2008) Improved segmental isotope labeling methods for the NMR study of multidomain or large proteins: application to the RRM of Npl3p and hnRNP L. *J. Mol. Biol.* 375, 151-164. (pTWIN1) PubMedID: [17936301](#).

Zhao, W. et al. (2008) An efficient on-column expressed protein ligation strategy: application to segmental triple labeling of human apolipoprotein E3. *Protein Sci.* 17, 736-747. (pTYB1) PubMed ID: [18305193](#).

Ottesen, J. et al. (2008) An amalgamation of solid phase peptide synthesis and ribosomal peptide synthesis. *Biopolymers*, 90, 406-414. (pTXB1) PubMed ID: [17636509](#).

2007

Hauser, P. S., et al. (2007) Expressed protein ligation using an N-terminal cysteine containing fragment generated *in vivo* from a pelB fusion protein. *Protein Expr. Purif.* 54, 227-233. (pTYB1) PubMed ID: [17493830](#).

Karow, A., et al. (2007) Authentic interdomain communication in an RNA helicase reconstituted by expressed protein ligation of two helicase domains. *FEBS J.* 2, 463-473. (pTWIN1) PubMedID: [17229151](#).

Kochinyan, S., et al. (2007) Use of intein-mediated phosphoprotein arrays to study substrate specificity of protein phosphatases. *Biotechniques*, 42, 63-69. (pTWIN2) PubMed ID: [17269486](#).

Paulick, M. G., et al. (2007) Synthetic analogues of glycosylphosphatidylinositol-anchored proteins and their behavior in supported lipid bilayers. *J. Am. Chem. Soc.* 129, 11543-11550. (pTXB1) PubMed ID: [17715922](#).

Sun, L., et al. (2007) Design, preparation and use of ligated phosphoproteins: A novel approach to study protein phosphatases by dot blot array, ELISA and Western blot assays. *Methods*, 42, 220-226. (pTXB1, pTWIN1) PubMed ID: [17532508](#).

2006

Kalia, J., et al. (2006) Reactivity of intein thioesters: appending a functional group to a protein. *Chembiochem.* 7, 1375-1383. (pTYB1) PubMed ID: [16897799](#).

Grant, J., et al. (2006) The N terminus of GTP gamma S-activated transducin alpha-subunit interacts with the C terminus of the cGMP phosphodiesterase gammasubunit. *J. Biol. Chem.* 281, 6194-6202. (pTXB1) Pubmed ID: [16407279](#).

Gottlieb, D., et al. (2006) Intein-mediated *in vitro* synthesis of lipidated Ras proteins. *Chem. Commun. (Camb)*. 21, 260-262. (pTYB2) PubMed ID: [16391727](#).

Eckenroth, B., et al. (2006) Semisynthesis and characterization of mammalian thioredoxin reductase. *Biochem.* 45, 5158-70. (pTYB3) PubMed ID: [16618105](#).

Lee, N.Y., et al. (2006) Structure and dynamics of the epidermal growth factor receptor C-terminal phosphorylation domain. *Protein Sci.* 15, 1142-52. (pTYB2) PubMedID: [16597832](#).

2005

Girish, A., et al. (2005) Site-specific immobilization of proteins in a microarray using intein-mediated protein splicing. *Bioorg. Med. Chem. Lett.* 15, 2447-51. (pTWIN) PubMedID: [15863295](#).

2004

Durek, T., et al. (2004) Synthesis of fluorescently labeled mono- and diprenylated Rab7 GTPase. *J. Am. Chem. Soc.* 126, 16368-16378. (pTYB1).

Durek T., Goody R.S and Alexandrov K. (2004) In vitro semisynthesis and applications of C-terminally modified Rab proteins. *Methods Mol Biol.* 283, 233-244. (pTWIN).

Hemelaar, J. et al. (2004) Specific and covalent targeting of conjugating and deconjugating enzymes of ubiquitin-like proteins. *Mol. Cell Biol.* 24, 84-95. (pTYB1).

Lue, R.Y., et al. (2004) Versatile protein biotinylation strategies for potential high-throughput proteomics. *J Am. Chem. Soc.* 126,1055-1062. (pTYB1, pTYB2).

Pezza, J.A., Allen, K.N. and Tolan, D.R. (2004) Intein-mediated purification of a recombinantly expressed peptide. *Chem Commun (Camb)*. 21, 2412-2413. (pTWIN1, pTYB1).

Romanelli, A., et al. (2004) Semisynthesis of a segmental isotopically labeled protein splicing precursor: NMR evidence for an unusual peptide bond at the N-extein-intein junction. *Proc. Natl. Acad. Sci. USA.* 101, 6397-6402. (pTXB1).

Takeda, S., et al. (2004) Site-specific conjugation of oligonucleotides to the C-terminus of recombinant protein by expressed protein ligation. *Bioorg. Med. Chem. Lett.* 14(10):2407-10. PubMedID: [15109622](#) (pTWIN1).

2003

David, R., et al. (2003) Semisynthesis and application of carboxyfluorescein-labelled biologically active human interleukin-8. *Biol. Chem.* 384, 1619-1630. (pTXB1)
PubMed ID: [14719805](#).

Hong, S.H., and Maret, W. (2003) A fluorescence resonance energy transfer sensor for the beta-domain of metallothionein. *Proc. Natl. Acad. Sci. USA.* 100, 2255-2260. (pTYB11).

Scheibner, K.A., Zhang, Z. and Cole, P.A. (2003) Merging fluorescence resonance energy transfer and expressed protein ligation to analyze protein-protein interactions. *Anal. Biochem.* 317, 226-232. (pTYB2).

Sinars, C.R., et al. (2003) Structure of the large FK506-binding protein FKBP51, an Hsp90-binding protein and a component of steroid receptor complexes. *Proc. Natl. Acad. Sci. USA.* 100, 868-873. (pCYB4).

Yee, C.S., et al. (2003) Generation of the R2 subunit of ribonucleotide reductase by intein chemistry: insertion of 3nitrotyrosine at residue 356 as a probe of the radical initiation process. *Biochemistry*, 42, 14541-14552. (pTYB3).

Zhang, Z., et al. (2003) The role of C-terminal tyrosine phosphorylation in the regulation of SHP-1 explored via expressed protein ligation. *J. Biol. Chem.* 278, 4668-4674. (pTYB2).

Zuberek, J., et al. (2003) Phosphorylation of eIF4E attenuates its interaction with mRNA 5' cap analogs by electrostatic repulsion: intein-mediated protein ligation strategy to obtain phosphorylated protein. *RNA*, 9, 52-61. (pTXB1).

2000

Cotton, G.J. and Muir, T.W. (2000). Generation of a dual-labeled fluorescence biosensor for Crk-II phosphorylation using solid-phase expressed protein ligation. *Chem. Biol*, 4, 253-261.

Structure studies

2003

Rak, A., et al. (2003) Structure of Rab GDP-dissociation inhibitor in complex with prenylated YPT1 GTPase. *Science*, 302, 646-50. (pTWIN1)

Walters, K.J., et al. (2003) DNA -repair protein hHR23a alters its protein structure upon binding proteasomal subunit S5a. *Proc. Natl. Acad. Sci. USA*. 100, 12694-12699. (pTWIN)

2002

Camarero, J.A., et al. (2002) Autoregulation of a bacterial sigma factor explored by using segmental isotopic labeling and NMR. *Proc. Natl. Acad. Sci. USA*. 99, 8536-8541. (pTXB1)

2001

Wu, J.W., et al. (2001) Crystal structure of a phosphorylated Smad2. Recognition of phosphoserine by the MH2 domain and insights on Smad function in TGF-beta signaling. *Mol Cell*. 8, 1277-1289. (pTXB1)

For additional references see [InBase](#), the Intein registry website.